Operation Manual

Shodex KW404-4F

J-24

(Please read this manual carefully before using the column to keep its good performance and life.)

1. Introduction

Shodex KW404-4F is a size exclusion chromatography column developed for separation of hydrophilic substances with low molecular weight, such as small proteins and peptides.

2. Instructions in handling <Important>

Caution!

* Take notice of keeping instructions about the solvents and the reagents used with the column not to occur problems related to losing your health or leaking.

Attention!

* Use the column within the regular range of flow rate, pressure and temperature. There is a danger of deteriorating the performance when it is handled beyond the permissible range even for a short time. See the clause "Usable conditions" about the permissible range.

3. Specifications

Column size	: 4.6 mm I.D. × 300 mm L.
Column material (inside)	: SUS 316
Packing material	: Porous silica particle (3 μm) bound with hydrophilic polymer
In-column solvent (initial)	: Water
Number of theoretical plates	: ≥25,000 per column

4. Usable conditions

Flow rate	: ≤0.5 ml/min
	* It should be 0.15 ml/min or less for replacing the in-column solvent.
Pressure	: ≤10MPa per column
Temperature	: 5~45 °C
pH	: 3.0~7.5
Eluent	: An aqueous solution of salt or a buffer solution is usually used.
	<Representative salts> sodium chloride, sodium nitrate, sodium sulfate, potassium sulfate, ammonium sulfate.
	<Representative buffers> phosphate buffer, tris buffer, acetate buffer, citrate buffer.
	* Methanol or acetonitrile can be added up to 100 vol%.

Attention!

- 1) Do not remove the end fittings of the column under any circumstances.
- 2) Do not make a strong impact on the column: such as hitting or dropping on the floor.
- 3) Replace the solvent in the chromatograph with the eluent to be used before connecting the column.
- 4) Connect the column so that the flow direction corresponds to the arrow mark on the tag.
- 5) When the column is not used for a month or more, replace the in-column solvent with water, close each end with a stopper, and store it at room temperature.
- 6) Adjust the total salt concentration of eluent to not more than 0.5 M.
- 7) Filtrate the sample with a disposable filter (0.45 μm) to prevent deterioration by adsorbing insoluble matters.
- 8) When methanol or acetonitrile is added to an aqueous solution of salt or a buffer solution, pay attention not to separate any salt after mixing.

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