Shodex[®]

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1. Introduction

Shodex SUGAR KS-800 series are packed with ionexchange resin gels of sulfonated styrene-divinylbenzene copolymer. They are Designed for High-speed and high-performance liquid chromatography. They are capable of separating carbohydrates such as monosaccharides, disaccharides, oligosaccharides and sugaralcohols, as well as nonionic water-soluble substance, with distilled water or 0.001-Mol sodium hydroxide aqueous solution used as the mobile phase. The columns demonstrate an excellent performance, particularly, in high-speed analysis of saccharides present in living bodies, food, pharmaceuticals and industrial wastewater.

2. Specifications

1) Nomenclature, exclusion limit and number of theoretical plates

Nomenclature	Exclusion limit	Number of theoretical plate /30cm ¹⁾
SUGAR KS-801	1×10^3	17,000 minimum
SUGAR KS-802	1×10^4	17,000 minimum
SUGAR KS-803	5×10^4	17,000 minimum
SUGAR KS-804	4×10^5	17,000 minimum
SUGAR KS-805	5×10^6	9,000 minimum
SUGAR KS-806	5×10^7 (est)	9,000 minimum

NOTE 1) : Specimen, Ethylene glycol(KS-801~803), Glucose(KS-804~806); mobile phase distilled water; flow rate 1mL/min; column temperature 50°C.

- 2) Size: ID 8mm, length 300mm
(exclusive of KS-G of 6mm in ID and 50mm in length.)
- 3) End fitting: Internally-threaded type, No.10-32 UNF.
- 4) Material: SUS 316
- 5) Packing: Strong cation-exchange resin gels.
- 6) In-column mobile phase: Distilled water.

4 Heat the column to 40 – 50°C and start the pump.

Do not let the air into the column during the connection.

Cautions

arrow mark on the column will point to the direction of the flow.

3) Connect the column to the chromatograph in such a way that the

2) Set the flow rate at 1mL/min.

replace the water-soluble solvent with the mobile phase.

water-soluble organic solvent, such as acetone or ethanol, and then

In replacing a water-insoluble organic solvent, replace it first with a

Cautions

the solvent in the chromatograph with the mobile phase.

1) Before mounting the column on the liquid chromatograph, replace

4. Mounting and start-up

the adsorptivity.

ethanol in an amount of 20% maximum to the mobile phase to reduce

specimens are too slow to flow out of the column. In such a case, add

Adsorbed sometimes by the packing, molecules of highly hydrophobic

Note

Use of solvent degassing device will facilitate the degassing work.

pressure being simultaneously reduce by an aspirator.

shaking in an ultrasonic washing vessel with the ambient air

degas it by, for instance, heating it about 60°C in a hot water bath and

by passing it through a 0.45 μm filter into a bottle and thoroughly

2) Remove extraneous substances and insolubles from the mobile phase

as the mobile phase.

1) Use distilled water or 0.001-Mol sodium hydroxide aqueous solution

3. Mobile phase

9) Maximum working temperature: 85°C.

1.5mL/min

8) Maximum working flow rate when column is heated to 50 – 85°C:

7) Maximum working pressure per column: 5MPa

6 . Dismounting and storage
 1) Set the flow rate at 0.2mL/min, turn off the heater and let the column stand as is until it cools down to room temperature.

4) When the specimen contains proteins, add 2 or 3 drops of distilled water in which sulfosalicylic acid is dissolved at the rate of 20g per deciliter, to 5mL of the specimen, stir and pass it through a 0.45 μm filter to remove the proteins.
 3) Pass the specimen through a 0.45 μm filter to remove extraneous substances and insolubles. Use of the disposable filter unit Shodex DT is recommended.
 2) When the specimen is liquid and contains an organic solvent such as alcohol, dilute it with the mobile phase so as to reduce concentration of the organic solvent to 20% maximum; otherwise, the column performance will deteriorate.
 Do not use any other solvent; otherwise, the solute sometimes settles in the mobile phase to plug the column.

5 . Pre-treatment of specimen
 1) When the specimen is solid, dissolve it, using some of the degassed mobile phase.
 Note

Working temperature	50 - 85°C	1.5mL/min
Maximum working flow rate	30 - 49°C	1.0mL/min
	29°C or below	0.5mL/min
	!!!) Do not abruptly change the column pressure or flow rate.	
	iv) Do not heat the column higher than 85°C.	

temperature, as follows:
 !!) The Maximum working flow rate depends on the working the column will be ruined.

!) Keep the column pressure below 5MPa per column; otherwise,
Caution

- 2) Stop the pump and dismount the column.
- 3) Cap both ends of the column, place it back in the same case in which it was delivered from the manufacturer, and store the case in a place where the temperature does not vary more than $\pm 5^{\circ}\text{C}$ from the ambient temperature.

Caution

!) Do not remove the end fittings of the column; otherwise, its performance will deteriorate or the column may be rendered unserviceable.

!!) Do not let the in-column solvent freeze during storage; otherwise, the column will swell and become unserviceable.

!!!) Storage at abnormally high temperatures will also cause the column to swell and become unserviceable.

7. Band broadening

Adsorption of organics and metal ions by the packing sometimes causes the bands of fructose and glucose to be broadened.

If such adsorption occurs, pass 0.2N NaOH aqueous solution through the column at a rate of 0.5mL/min to desorb them and completely replace the solution with distilled water. If, however, the band broadening is moderate, injection of 40 μL of 1N NaOH aqueous solution into the column will suffice to carry out the desorption.

8. Guard column and column jacket

Use of precolumn, SUGAR KS-G, is recommended immediately before the column to protect the packing from contamination with pollutants or readily adsorbable substances, such as metal ions and amino compounds.

9. Performance test

The performance of the column can be checked by calculating the number of theoretical plates under the following conditions.

- 1) Specimen and Injection amount

KS-801	1.0% ethylene glycol aqueous solution	5 μ L
KS-802	2.5% ethylene glycol aqueous solution	5 μ L
KS-803	0.5% ethylene glycol aqueous solution	15 μ L
KS-804	0.8% glucose aqueous solution	5 μ L
KS-805	0.4% glucose aqueous solution	5 μ L
KS-806	0.8% glucose aqueous solution	5 μ L
- 2) Mobile phase: Distilled water
- 3) Flow rate: 1.0mL/min
- 4) Column temperature: 50°C
- 5) Detector: High-sensitivity refractive index detector
such as Shodex RI

6) Calculation formula

$$NTP=5.54 \times (Rt/W)^2$$

where NTP: Number of theoretical plate
Rt: Retention time
W: Peak half width