



**Shodex™**  
CAPTURE THE ESSENCE

# LC/UV/MS analysis of monitoring bioethanol manufacturing process with using polymer-based multi-solvent SEC column

Junji Sasuga<sup>2</sup>, Melissa Turcotte<sup>1</sup>, and Ronald Benson<sup>1</sup>

<sup>1</sup>Showa Denko America, Inc., 420 Lexington Avenue, Suite 2335A, New York, NY 10170, USA

<sup>2</sup>Showa Denko K.K., 5-1 Ohgimachi Kawasaki-ku, Kawasaki, 210-0867, Japan

Contact us at [support@shodex.net](mailto:support@shodex.net).  
For more information, visit [www.shodex.net](http://www.shodex.net)

## Abstract

The production of bioethanol begins with pyrolysis of bagasse, enzymatic conversion into saccharides, and finally fermentation. In order to improve productivity, saccharides as well as inhibitors (e.g. organic acids and furfurals) of fermentation should be measured.


Typically, saccharides are analyzed by ligand exchange chromatography, while organic acids and furfurals are analyzed separately by ion exclusion/reversed-phase chromatography. In order to simplify the analysis, a simultaneous LC/UV/MS method was studied.

A multi-solvent size exclusion chromatography (SEC) column with poly vinyl alcohol base packing material, Shodex Asahipak GF-210 HQ was used. A mixed standard solution containing glucose, cellooligosaccharides (2-5 saccharide units) and various organic acids was analyzed. The saccharides eluted from the column in accordance with SEC mode, followed by elution of organic acids and furfurals. Retention of the organic acids was attributed to hydrophobic interaction with the packing material rather than SEC mode. Furfurals were detected more sensitively by UV than by MS.

A sample obtained by pyrolysis of bagasse and enzymatic conversion to saccharides was also analyzed. The analysis showed the standard substances as well as xylooligosaccharides, organic acids originated from hemicellulose, and benzaldehydes originated from lignin.

Our method using simultaneous LC/UV/MS is effective for analysis of bioethanol production.

## Manufacturing process of bioethanol



- The components (organic acids, furfurals, and so on) produced in thermal hydrolysis interfere with fermentation, so these products should be analyzed in order to control them.
- Saccharides and these products have been ordinarily measured at different conditions (SEC or ion exclusion/reversed-phase chromatography).

Simultaneous analysis of all compounds in one condition is very efficient.

## Aqueous and organic SEC columns Asahipak GF-HQ series


### Product specifications and line up

| Product name | Exclusion Limit (Pululan) | Plate Number (TP/Column) | Particle Size (μm) | Column Size (I.D. x Length (mm)) |
|--------------|---------------------------|--------------------------|--------------------|----------------------------------|
| GF-210 HQ    | 9,000                     | ≥ 19,000                 | 5                  | 7.5 x 300                        |
| GF-310 HQ    | 40,000                    | ≥ 19,000                 | 5                  | 7.5 x 300                        |
| GF-510 HQ    | 300,000                   | ≥ 19,000                 | 5                  | 7.5 x 300                        |
| GF-710 HQ    | (10,000,000)*             | ≥ 11,000                 | 9                  | 7.5 x 300                        |
| GF-7M HQ     | (10,000,000)*             | ≥ 13,000                 | 9                  | 7.5 x 300                        |
| GF-1G 7B     | -                         | (guard column)           | 9                  | 7.5 x 50                         |


Packing material : Poly vinyl alcohol particle **Polymer based**  
Housing : Stainless steel  
Usable temp. : 4~60 °C (recommend 25~40 °C)  
Usable pH range : 2~9  
Max flow rate : 1.0 mL/min (recommend 0.4~0.6 mL/min)  
Max pressure : 7 MPa (GF-210 HQ and GF-310 HQ)  
Usable solvents : Water, Methanol, Ethanol, Acetonitrile, THF, Acetone, DMF, Chloroform, and DMSO (0~50%)

Aqueous and various organic solvents available

### Calibration curves




### Separation of maltooligosaccharides




- GF-210 HQ is suitable for the separation of oligosaccharides than GF-310 HQ.
- Hydrophobic interaction with the packing material is available for the separation.

### Utilizing of hydrophobic interaction



## [Result 1] Simultaneous analysis of the mixed standard


### Mass chromatograms of the mixed standard




Sample : 5 μg/mL each (in water), 20 μL  
Instrument : Shimadzu Nexera / LCMS-8030  
Column : shown above each data  
Eluent : (A) Water/Formic acid = 100/0.1 (v/v)  
(B) Acetonitrile  
Isocratic ; (B%) 20%  
Flow rate : 0.6 mL/min  
Column temp. : 50 °C  
Detector : PDA & ESI-MS (SIM)

| Peak No. | Compound      | Monitored m/z |
|----------|---------------|---------------|
| 1        | Cellopentaose | 827(-)        |
| 2        | Cellobiose    | 665(-)        |
| 3        | Cellobiose    | 503(-)        |
| 4        | Glucose       | 341(-)        |
| 5        | Xylobiose     | 281(-)        |
| 6        | Xylose        | 149(-)        |
| 7        | Glyceric acid | 105(-)        |
| 8        | Glycolic acid | 75(-)         |
| 9        | Lactic acid   | 89(-)         |
| 10       | Malic acid    | 133(-)        |
| 11       | Succinic acid | 117(-)        |
| 12       | Tartaric acid | 149(-)        |
| 13       | Adipic acid   | 145(-)        |
| 14       | Citric acid   | 191(-)        |
| 15       | Malonic acid  | 103(-)        |
| 16       | Fumaric acid  | 115(-)        |
| 17       | Pyruvic acid  | 87(-)         |
| 18       | Maleic acid   | 115(+)        |
| 19       | 5-HMF         | 127(+)        |
| 20       | Furfural      | 97(+)         |

### UV chromatogram of the mixed standard



### Calibration curves




High linearity & high sensitivity (0.05 μg/mL detectable-each)

- Using Asahipak GF-210 HQ, saccharides were separated by SEC mode, and organic acids and furfurals were separated by reverse-phase mode after the elution of saccharides.
- But saccharides and organic acids were eluted in the same zone, using two analysis columns for organic acids from other manufacturers.
- Furfurals were sensitively detected by UV.

• LC/UV/MS using GF-210 HQ is very useful for simultaneously analyzing saccharides and other compounds such as organic acids and furfurals.

## [Result 2] Analysis of the actual sample

### Chromatograms of the saccharification sample




• Three kinds of benzaldehydes known as inhibitors of fermentation derived from lignin were detected by UV and MS.

• The main two unknown peak were detected at 58min and 79min on UV Chromatogram. These components may be novel benzaldehydes.

• This method was applicable for analysis of the thermal hydrolysis sample and saccharification sample (The data of the thermal hydrolysis sample is not shown here).

• Analysis time is estimated to become shorter by increasing the ratio of acetonitrile.



## Conclusion

LC/UV/MS method with GF-210 HQ enables simultaneously analyzing saccharides, organic acids, furfurals, and benzaldehydes.

This method provides an efficient and sensitive analysis of the samples of bioethanol manufacturing process.