

LC/MS Method Development for Various Oligonucleotides Using Polymer-Based HILIC Columns Having Diol Group



Shodex™
CAPTURE THE ESSENCE

Eiji Kagawa², Hisako Sakuma², Junji Sasuga², Kanna Ito², Leah Sullivan¹, and Ronald Benson¹

¹Showa Denko America, Inc., 420 Lexington Avenue, Suite 2335A, New York, NY 10170, USA

²Showa Denko K.K., 5-1 Ohgimachi Kawasaki-ku, Kawasaki, 210-0867, Japan

Introduction

Development and quality control of oligonucleotide therapeutics often require highly selective and sensitive analytical methods. An LC/MS method is one of them. The most often used LC method for the oligonucleotide analyses have been reversed phase with a use of ion-pair reagents. However, the ion-pair reagents have a tendency to remain on the LC system and cause problems.

In this work, we used Shodex™ HILICpak™ VN-50 series columns, polymer-based HILIC columns modified with diol functional groups. The developed method enabled the LC/MS analysis of oligonucleotides without a need of an ion-pair reagent. Gradient condition was optimized with a mixture of acetonitrile and relatively low-concentration volatile basic solvent.

Shodex™ HILICpak™ VN-50 Series Columns

Product name	Column size (mm)
	I.D. x Length
HILICpak VN-50 4D	4.6 x 150
HILICpak VN-50G 4A	4.6 x 10
HILICpak VN-50 2D	2.0 x 150
HILICpak VN-50G 2A	2.0 x 10
HILICpak VN-50 1D	1.0 x 150


※VN-50 1D will be available soon

< Features >

- Suitable for oligosaccharides and oligonucleotides analysis using HILIC mode
- Retains hydrophilic compounds without derivatization nor ion-pair reagent

< Specifications >

- Base Material: Spherical porous particles of polyvinyl alcohol modified with diol functional group
- Particle Size: 5 μm
- Pore Size: 100 Å
- Column Housing: PEEK
- Usable pH: 2 - 13 Workable in alkaline condition
- Usable temperature: 4 - 60 °C
- Applicable solvents: Mixture of water and acetonitrile or water and methanol of any ratio. Up to 100 % water, acetonitrile, or methanol.



Method

< Analytical Condition >

Instrument : Shimadzu Nexera / LCMS-8030 Plus
 Column : Shodex HILICpak VN-50 2D (2.0 mm I.D. x 150 mm)
 Shodex HILICpak VN-50 1D (1.0 mm I.D. x 150 mm)
 Eluent : (A) 50mM HCOONH₄ aq.(pH 9.8) / (B) CH₃CN
 Linear gradient (High pressure)
 Flow rate : 0.3 mL/min (VN-50 2D), 0.1 mL/min (VN-50 1D)
 Detector : PDA (190 - 350 nm) & ESI-MS SIM(-)
 Column temp.: 40 °C
 Injection vol. : 1 μL


No ion-pairing reagent nor HFIP used

< Sample >

- Synthetic phosphorothioated oligo-RNA (Cartridge purified)
1. 19 mer U*A*C*C*G*A*U*U*A*A*G*C*G*A*A*G*U*U*U
 2. 20 mer A*U*A*C*C*G*A*U*U*A*A*G*C*G*A*A*G*U*U*U
 3. 21 mer G*A*U*A*C*C*G*A*U*U*A*A*G*C*G*A*A*G*U*U*U
- Synthetic phosphorothioated oligo-DNA (salt-free)
4. 20 mer A*T*A*C*C*G*A*T*T*A*A*G*C*G*A*A*G*T*T*T
- Synthetic oligo-DNA (salt-free)
5. 20 mer ATACCGATTAAGCGAAGTTT


< Phosphorothioated Oligo-RNA Analysis >

Elute in order of smaller to larger oligo-RNAs by HILIC mode




< Repeatability >

Retention time deviations for all peaks: ≤ 0.5% RSD (n = 6)




Results and Discussion

< Comparisons of Oligonucleotides >




- Less retention for phosphorothioated oligonucleotides (higher hydrophobicity)
- Stronger retention for phosphorothioated oligo RNAs than phosphorothioated oligo-DNA (RNA with OH group at the ribose 2' position makes it more hydrophilic than DNA)




< Comparisons of Column I.D. Size >

Superior ESI-MS sensitivity with VN-50 1D



< Calibration Curves >

- Good overall linearity
- Need further work in the ng/mL level quantification



Conclusions

The Shodex™ HILICpak™ VN-50 series, polymer-based diol type HILIC columns, demonstrated highly selective LC/UV/ESI-MS measurements. The oligonucleotides were separated and eluted in the order of smaller to higher degree of polymerization. The gradient setting with a mixture of 50-mM aqueous ammonium formate solution and acetonitrile enabled the analysis without a need of ion-pair reagents for reversed phase mode or highly concentrated salt for ion exchange mode generally used. The developed method will also be beneficial for preparative purposes, as it does not require de-salting process. Furthermore, 1.0 mm I.D. column showed an improved sensitivity for the synthetic phosphorothioated oligonucleotides MS measurements. Because of their stabilities compared to regular oligonucleotides, phosphorothioated oligonucleotides are more often chosen for oligonucleotide therapeutics. The results presented in this work proved that the VN-50 series columns are expected to be effective tools for the development and quality control of oligonucleotide therapeutics.